



CiFLAB
SOLUTIONS

Ascent™

Features and Benefits:

The Ascent™ frames are designed for use as a single table or a back-to-back station and features rounded uprights that are used for both single or double sided applications.

Adjustment is accomplished by removing bolts holding the top frame to the uprights and repositioning. Machine screws are removed in front legs and repositioned.

Shelves are available in plastic laminate, painted metal, phenolic resin*, epoxy resin, and wood and adjustable in 1-inch increments

Optional shelf retainer lips are available in metal cut-out, single rail and triple rail

Accessories include apparatus bars, pneumatic height-adjustable flat screen monitor arms* and keyboard trays*

Horizontal utility manifolds meet electrical requirements*

Mobile casework options available in plastic laminate, painted metal* or wood

Manual* height-adjustable work surface

Full-height utility uprights act as a service chase for air, vacuum and gas

*Option shown on drawing

