



C-Frame System

Product Description and Characteristics:

C-Frames are affordable solutions, built to last. Our C-Frame systems are ordered as starters or adders in three heights: bench (nominal 36-inches), ADA (nominal 34-inches) and desk (nominal 30-inches). Every starter and adder unit can be designed with open backs or with 1/2-inch removable thermally fused closure panels available in six standard colors. When using with under counter refrigerators, contact CiF Lab Solutions for custom sizes.

Load Rating:

The total load rating is 600 LBS. and based on evenly distributed loads.

Finish:

All steel frame components shall be finished using a baked-on powder-coat finish. Standard colors include black, grey, light neutral, and white.

Fabrication:

C-Frame Systems are shipped knocked down on pallets with all accessory items. Suspended and mobile casework items are shipped assembled, but need to be installed on the C-Frame at final destination.

Warranty:

All standard products manufactured by CiF Lab Solutions carry a two (2) year guarantee and Limited Warranty from the date of substantial completion to the original owner against defective material and workmanship.

