



CiFLAB
SOLUTIONS


Carts - Manual Height Adjustable or Fixed Height

Features and Benefits:

The fixed height Cart is an economical solution to mobility. Carts are available with leveling glides or swivel casters.

Shelf brackets are available as angle*
or book-end style

Shelves are available in plastic laminate,
painted metal*, phenolic resin, epoxy
resin, and wood and infinitely adjustable



Optional Manual Height Adjustable Work
surface

Optional shelf retainer lips are available
in metal cut-out, single rail and
triple rail*

Horizontal gas and electric manifolds*

Available with leveler glides* and
swivel casters*

Work surface material options include
phenolic resin, thermally fused,
metal, epoxy resin, wood, and plastic
laminate*

*Option shown on drawing