




CiFLAB
SOLUTIONS

S-Line Series

Our durable metal casework line is designed to meet the most demanding of specifications. Robust construction, quality hardware and exceptional fit and finish make our S-line Series of metal casework a true value. We offer many features, such as full extension slides, removable backs and removable bottoms as standard construction. Our metal casework is available in inset overlay, flush overlay and flush overlay with wood fronts.



Door and drawer front construction allows for Inset Steel, Flush Overlay Steel or Flush Overlay Wood (on the steel cabinet body)



Full extension, ball bearing, plated drawer slides

Stainless steel, 5 knuckle institutional type hinge

Door and drawer heads are "double pan" construction with sound deadening interior

Standard pulls are flush extruded aluminum (with capability to be a label holder) or an optional aluminum bar style pull is available

Full depth adjustable shelf that is adjustable on 1" centers

Adjustable doors with full access opening

Removable backs

Removable bottom panned up to contain spills



SEFA (Scientific Equipment Furniture Association) . . . as a SEFA member, CiF builds their standard cabinets in compliance with all SEFA requirements.