



CiFLAB  
SOLUTIONS

## Tech Center

### Features and Benefits:

The Tech Center is available as a fixed height workstation or with optional leg extensions that allow for 8-inches of height adjustment. The leg extensions are inserted inside the workstation legs and secured with machine screws. For increased mobility, order your Tech Center with swivel casters.

Tech Center uprights are deck-mounted and available as aluminum infinitely adjustable (36- or 42-inch height uprights) and steel slotted adjustable (36- or 48-inch height uprights) in nominal 1-inch increments. When ordering workstations with a width over 48-inches, order the uprights that include the center upright.

Because you may have special applications and requirements, we offer several optional accessories, including gas, air and vacuum manifolds, flat-screen monitor arms, CPU holders and more.

---

Shelf brackets are available as angle or book-end\* style

---

Shelves are available in plastic laminate, painted metal\*, phenolic resin, epoxy resin, and wood and adjustable in 1-inch increments

---

Optional shelf retainer lips are available in metal cut-out\*, single rail, and triple rail

---

Available with leveler glides, swivel casters or leg extensions with leveler glides\*

---

\*Option shown on drawing

