



Contents

Metal Base Cabinet Details (S-Line Series).....	8.02
Metal Base Cabinets.....	8.03-8.10
Metal Sink Base Cabinets	8.11-8.12
Metal Hood Base Cabinets	8.13
Metal Acid Storage Base Cabinets.....	8.14
Metal Flammable Storage Base Cabinets.....	8.15



CiF LAB
SOLUTIONS


S-Line Series

Our durable metal casework line is designed to meet the most demanding of specifications. Robust construction, quality hardware and exceptional fit and finish make our S-line Series of metal casework a true value. We offer many features, such as full extension slides, removable backs and removable bottoms as standard construction. Our metal casework is available in inset overlay, flush overlay and flush overlay with wood fronts.



SEFA (Scientific Equipment Furniture Association)... as a SEFA member, CiF builds our standard cabinets in compliance with all SEFA requirements.

General Storage Base Cabinets



CAT. NO

B0000

DESCRIPTION

Open Cabinet
• One shelf

HEIGHT: INCHES


29 (Sitting) 33 (ADA) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24 30 36 42 48



CAT. NO

B0110 (Hinged L) B0120 (Hinged R)

DESCRIPTION

Storage w/Door
• One shelf

HEIGHT: INCHES


29 (Sitting) 33 (ADA) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

12 15 18 24




B0100

Storage w/Doors
• One shelf

29 (Sitting) 33 (ADA) 35 (Standing)

22

30 36 42 48




B3000

Door/Drawer Storage
• One shelf on 35" H only

29 (Sitting) 33 (ADA) 35 (Standing)

22

30 36 42 48



CAT. NO

B3100

DESCRIPTION

Door/Drawer Storage
• One shelf on 35" H only

HEIGHT: INCHES


29 (Sitting) 33 (ADA) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48




B3110 (Hinged L) B3120 (Hinged R)

Door/Drawer Storage
• One shelf on 35" H only

29 (Sitting) 33 (ADA) 35 (Standing)

22

15 18 24



B3400


Door/Drawer Storage
• Double hinged doors and four half high drawers, one shelf

29 (Sitting) 35 (Standing)

22

30 36 42 48

General Storage Base Cabinets



CAT. NO

B3410 (Hinged L) B3420 (Hinged R)

DESCRIPTION

Door/Drawer Storage

- One hinged door and two drawers, one shelf on 35"H only

HEIGHT: INCHES


29 (Sitting) 33 (ADA) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24



B3430


Door/Drawer Storage

- Double hinged doors, two left half high drawers and one right drawer, one shelf

35 (Standing)

22

30 36 42 48



B3440


Door/Drawer Storage

- Double hinged doors, left drawer and two right half high drawers, one shelf

35 (Standing)

22

30 36 42 48



CAT. NO

B4030

DESCRIPTION

Three Drawer Storage

- Three drawers

HEIGHT: INCHES


29 (Sitting)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24



B4040

Four Drawer Storage

- Four drawers

33 (ADA) 35 (Standing)

22

15 18 24 30 36 42 48



B4070


Seven Drawer Storage

- Seven drawers

33 (ADA)

22

15 18 24 30 36 42 48



CAT. NO

B4080

DESCRIPTION

Eight Drawer Storage

- Eight drawers

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24 30 36 42 48



B4160


Six Drawer Storage

- Six drawers

29 (Sitting)

22

18 24 30 36 42 48



B4240

Five Drawer Storage


- Three drawers and two half-width drawers

35 (Standing)

22

30 36 42 48

General Storage Base Cabinets



CAT. NO

B4290

DESCRIPTION

Three Drawer Storage
• Three drawers

HEIGHT: INCHES


29 (Sitting)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24




B4330

Six Drawer Storage
• Six drawers (2 x 3)

29 (Sitting)

22

30 36 42 48




B4340

Eight Drawer Storage
• 8 drawers (2 x 4)

33 (ADA) 35 (Standing)

22

30 36 42 48



CAT. NO

B4370

DESCRIPTION

Fourteen Drawer Storage
• Fourteen drawers (2 x 7)

HEIGHT: INCHES


33 (ADA)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48




B4390

Four Drawer Storage
• Two drawers and two half high drawers

29 (Sitting)

22

15 18 24




B4420

Door / Drawer Storage
• Two full-width drawers, two left drawers, one full and two half high right drawers

35 (Standing)

22

30 36 42 48



CAT. NO

B4430

DESCRIPTION

Drawer Storage
• Three drawers and two half high drawers

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24 30 36 42 48



B4440


Drawer Storage
• Two drawers and four half high drawers

35 (Standing)

22

15 18 24 30 36 42 48

General Storage Base Cabinets



CAT. NO

B4460

DESCRIPTION

Drawer Storage

- Two full width drawers, two left drawers, and four right half high drawers

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B4470


Drawer Storage

- Four full width half high drawers, two left drawers and two right drawers

35 (Standing)

22

30 36 42 48



B4480


Drawer Storage

- Two full-width drawers and four drawers

35 (Standing)

22

30 36 42 48



CAT. NO

B4500

DESCRIPTION

File Drawer Storage

- Two file drawers

HEIGHT: INCHES

35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24



B4510

File Drawer Storage

- One file and one drawer

B4540


File Drawer Storage

- One file and two drawers

29 (Sitting) 33 (ADA) 35 (Standing)

22

15 18 24



CAT. NO

B4550

DESCRIPTION

File Drawer Storage

- One file, one drawer and two half high drawers

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24



B4610


Drawer Storage

- One full width drawer and six half width drawers

35 (Standing)

22

30 36 42 48



B4620

Drawer Storage


- Six half width drawers and one full width top drawer

35 (Standing)

22

30 36 42 48

General Storage Base Cabinets



CAT. NO

B4630

Drawer Storage

- One full width drawer and four half width drawers

DESCRIPTION

HEIGHT: INCHES


29 (Sitting)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B4640


Drawer Storage

- Four half width drawers and four half high drawers

29 (Sitting)

22

30 36 42 48



B4650


Drawer Storage

- Six half width drawers and four half high drawers

35 (Standing)

22

30 36 42 48



CAT. NO

B4660

Drawer Storage

- Four half width drawers and 8 half high drawers

DESCRIPTION

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B5410


Door/Drawer Storage

- Left hinged door/drawer and four right drawers, one shelf

35 (Standing)

22

30 36 42 48



B5420


Door/Drawer Storage

- Right hinged door/drawer and four left drawers, one shelf

35 (Standing)

22

30 36 42 48



CAT. NO

B6330

Door/Drawer Storage

- Left hinged door, full-width top drawer and three right drawers, one shelf

DESCRIPTION

HEIGHT: INCHES


33 (ADA) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B6340


Door/Drawer Storage

- Right hinged door, full-width top drawer and three left drawers, one shelf

33 (ADA) 35 (Standing)

22

30 36 42 48



B6370

Door/Drawer Storage


- Left hinged door, two half high drawers and four right drawers, one shelf

35 (Standing)

22

30 36 42 48

General Storage Base Cabinets



CAT. NO

B6380

DESCRIPTION

Door/Drawer Storage

- Right hinged door, two half high drawers and four left drawers, one shelf

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B6410


Door/Drawer Storage

- Left hinged door and four right drawers, one shelf

33 (ADA) 35 (Standing)

22

30 36 42 48



B6420


Door/Drawer Storage

- Right hinged door and four left drawers, one shelf

33 (ADA) 35 (Standing)

22

30 36 42 48



CAT. NO

B6450

DESCRIPTION

Door/Drawer Storage

- Left hinged door, two half high drawers and three right drawers

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B6460


Door/Drawer Storage

- Right hinged door, two left half high drawers, and three left drawers

35 (Standing)

22

30 36 42 48



B6630


Door/Drawer Storage

- Left hinged drawer, one left drawer, two right drawers and four half high drawers

35 (Standing)

22

30 36 42 48



CAT. NO

B6640

DESCRIPTION

Door/Drawer Storage

- Right hinged door, two left drawers, four half high drawers and one right drawer

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B6710


Door/Drawer Storage

- Left hinged door and seven right half high drawers

33 (ADA)

22

30 36 42 48



B6720

Door/Drawer Storage

- Right hinged door and seven left half high drawers

33 (ADA)


22

30 36 42 48



CiFLAB
SOLUTIONS

General Storage Base Cabinets



CAT. NO

B6880

Drawer Storage

- Four left drawers and seven right half high drawers

DESCRIPTION

HEIGHT: INCHES


35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B6890

Drawer Storage


- Four right drawers and seven left half high drawers

35 (Standing)

22

30 36 42 48

Corner Base Cabinets



CAT. NO

B0210

DESCRIPTION

Blind Corner

- Fixed panel on left and left hinged door on right

HEIGHT: INCHES


29 (Sitting) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48



B0220


Blind Corner

- Fixed panel on right and right hinged door on left

29 (Sitting) 35 (Standing)

22

30 36 42 48



B3210


Door/Drawer Blind Corner

- One shelf on 35" H only

29 (Sitting) 35 (Standing)

22

30 36 42 48



CAT. NO

B3220

DESCRIPTION

Door/Drawer Blind Corner

- One shelf on 35" H only

HEIGHT


29 (Sitting) 35 (Standing)

DEPTH

22

WIDTH: INCHES

30 36 42 48



CAT. NO

X1210 (Hinged L) X1220 (Hinged R)

DESCRIPTION

45 Degrees Corner Cabinet

- Angle corner
- Lazy Susan (LS000) ordered separately

HEIGHT: INCHES

29 (Sitting) 35 (Standing)


DEPTH: INCHES

22

WIDTH: INCHES

31

Sink and Utility Base Cabinets



CAT. NO

B2100

DESCRIPTION

Sink Cabinet w/Doors
• No Shelf

HEIGHT: INCHES


29 (Sitting) 33 (ADA) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

30 36 42 48




B2110 (Hinged L) B2120 (Hinged R)

Single Door Sink Cabinet
• No Shelf

29 (Sitting) 33 (ADA) 35 (Standing)

22

12 15 18 24



B2200

Sink Cabinet w/Doors
• No Shelf
• Louvered Sink Panel

29 (Sitting) 33 (ADA) 35 (Sitting)

22

30 36 42 48 52 58



CAT. NO

B2210 (Hinged L) B2220 (Hinged R)

DESCRIPTION

Single Door Sink Cabinet
• No Shelf
• Louvered Sink Panel

HEIGHT: INCHES


29 (Sitting) 33 (ADA) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24



B2750

Sink Cabinet Island End
• No Shelf, w/Kick Returns

29 (Sitting) 33 (ADA) 35 (Standing)

22

52 58




X1250 (Hinged L)

Corner Base Cabinet w/Sink Panel

29 (Sitting) 35 (Standing)

22

14



CAT. NO

X1260 (Hinged R)

DESCRIPTION

Corner Base Cabinet w/Sink Panel

HEIGHT: INCHES

29 (Sitting) 35 (Standing)

DEPTH: INCHES

22

WIDTH: INCHES

14




X1270 (Hinged L)

Corner Base Cabinet w/Louvered Sink Panel

29 (Sitting) 35 (Standing)

22

14



X1280 (Hinged R)

Corner Base Cabinet w/Louvered Sink Panel

29 (Sitting) 35 (Standing)


22

14



CiFLAB
SOLUTIONS

Sink and Utility Base Cabinets



CAT. NO

SS040

DESCRIPTION

Sink Support Assembly

- Designed for use with CiF Lab Solutions metal sink base cabinets only

HEIGHT: INCHES

DEPTH: INCHES




WIDTH: INCHES

12



CiF LAB
SOLUTIONS

Fume Hood Base Cabinets





CAT. NO	B2600	B2610 (Hinged L) B2620 (Hinged R)	B2650
DESCRIPTION	Fume Hood Cabinet • Double hinged doors with flush panel	Fume Hood Cabinet • Hinged door with flush panel	Fume Hood Cabinet • Double hinged door with flush panel
HEIGHT: INCHES	29 (Sitting) 35 (Standing)	29 (Sitting) 35 (Standing)	35 (Standing)
DEPTH: INCHES	22	22	22
WIDTH: INCHES	30 36 42 48 54	12 15 18 24	52 58



CiF LAB
SOLUTIONS

Acid Storage Base Cabinets



CAT. NO

SC230

SC231 (Hinged L)
SC232 (Hinged R)

DESCRIPTION

Acid Storage

- Double doors, no panel

Acid Storage

- One door, no panel
- Drawing is for concept only. Single door cabinet will be supplied

HEIGHT: INCHES

29* (Sitting) 35 (Standing)

29* (Sitting) 35 (Standing)

DEPTH: INCHES

18 22


18 22

WIDTH: INCHES

30 36 42 48 52 58

18 24

*Note: Spacer apron required for ADA installation.



CAT. NO

SCV00



DESCRIPTION

Vent Kit



CiFLAB
SOLUTIONS

Metal Flammable Storage Base Cabinets



CAT. NO

SC310

**SC311 (Hinged L)
SC312 (Hinged R)**

DESCRIPTION

**UL Listed Metal Flammable Storage Cabinet
(Manual Closing)**
• Double doors

**UL Listed Metal Flammable Storage Cabinet
(Manual Closing)**
• Single door

HEIGHT: INCHES

29* (Sitting) 35 (Standing)

29 (Sitting) 35 (Standing)

DEPTH: INCHES

18 22



18 22

WIDTH: INCHES

30 36 42 48 52

24

*Note: Spacer apron required for ADA installation.



CAT. NO

SC330

**SC331 (Hinged L)
SC332 (Hinged R)**

DESCRIPTION

**Non-UL Flammable Storage Base Cabinet
(Self Closing/Latching)**

**Non-UL Flammable Storage Base Cabinet
(Self Closing/Latching)**

HEIGHT: INCHES

29 35

29 35

DEPTH: INCHES

18

18

WIDTH: INCHES

30 36 42 48 52 58

30 36 42 48 52 58




CiFLAB
SOLUTIONS

Contents

Metal Wall Cabinets.....	9.02-9.04
--------------------------	-----------

General Storage Wall Cabinets



CAT. NO

W0000

DESCRIPTION

Open Storage

HEIGHT: INCHES


24 31 36 48

DEPTH: INCHES

13 16

WIDTH: INCHES

18 24 30 36 42 48




W0110 (Hinged L) W0120 (Hinged R)

Storage w/Door

24 31 36 48

13 16

18 24




W0100

Storage w/Doors

24 31 36 48

13 16

30 36 42 48



CAT. NO

W0840

DESCRIPTION

Storage w/Sliding Doors

HEIGHT: INCHES


24 31 36 48

DEPTH: INCHES

13 16

WIDTH: INCHES

30 36 42 48




W0310 (Hinged L) W0320 (Hinged R)

Storage w/Framed Glazed Door

24 31 36 48

13 16

18 24




W0300

Storage w/Hinged Framed Glazed Doors

24 31 36 48

13 16

30 36 42 48



CAT. NO

W0870

DESCRIPTION

Storage w/Sliding Framed Glazed Door

HEIGHT: INCHES


24 31 36 48

DEPTH: INCHES

13 16

WIDTH: INCHES

30 36 42 48



W0860

Storage w/Sliding Doors

24 31 36 48

13 16




30 36 42 48

Note:
Shelf quantity varies with height. Quantities below are typical for all wall metal cabinets.
24"H = 1 shelf
31"H = 2 shelves
36"H = 2 shelves
48"H = 3 shelves



CiFLAB
SOLUTIONS

General Storage Pass Through Wall Cabinets





CAT. NO	Y1200	W4000	W4870
DESCRIPTION	Corner Wall Case w/Open Angle Corner	Storage w/Sliding Doors	Storage w/Unframed Sliding Glass Doors
HEIGHT: INCHES	31 36	31	31 36
DEPTH: INCHES	13	13	13
WIDTH: INCHES	18	30 36 48	36 48



CiF LAB
SOLUTIONS

Corner Wall Cabinets



CAT. NO

Y1210 (Hinged L) Y1220 (Hinged R)

Y1270 (Hinged L) Y1280 (Hinged R)

DESCRIPTION

Angled Corner w/Door

Angled Corner w/Framed Glass

HEIGHT: INCHES

24 31 36 48

24 31 36 48

DEPTH: INCHES

13 16

13 16

WIDTH: INCHES

25 11/16 28 11/16

25 11/16 28 11/16




CiFLAB
SOLUTIONS

Contents

Metal Tall Cabinets	10.02-10.03
----------------------------------	--------------------

General Storage Tall Cabinets



CAT. NO

T0000

DESCRIPTION

Open Storage
• Five shelves

HEIGHT: INCHES


84

DEPTH: INCHES

16 21-1/4

WIDTH: INCHES

18 24 30 36 42 48



T0110 (Hinged L) T0120 (Hinged R)


Storage w/Door

• Hinged Single solid door, five shelves

84

16 21-1/4

18 24



T0100


Storage w/Doors

• Hinged Double solid doors, five shelves

84

16 21-1/4

30 36 42 48



CAT. NO

T0840

DESCRIPTION

Storage w/Sliding Doors
• Sliding solid doors, five shelves

HEIGHT: INCHES


84

DEPTH: INCHES

16 21-1/4

WIDTH: INCHES

30 36 42 48



T0310 (Hinged L) T0320 (Hinged R)


Storage w/Framed Glazed Door

• Hinged glazed door, five shelves

84

16 21-1/4

18 24



T0300

Storage w/Framed Glazed Doors

• Hinged Double glazed doors, five shelves

84


16 21-1/4

30 36 42 48



CiFLAB
SOLUTIONS

General Storage Tall Cabinets



CAT. NO

T0860

DESCRIPTION

Sliding Framed Glazed Doors
• Five shelves

HEIGHT: INCHES


84

DEPTH: INCHES

16 21-1/4

WIDTH: INCHES

30 36 42 48



T0840

Storage w/Sliding Doors

- Four adjustable shelves
- One centered fixed shelf
- Two sliding doors on overhead track

84

16 21-1/4


30 36 42 48



Contents

Metal Base Fillers, Panels & Accessories.....	11.02
Metal Wall Fillers	11.03
Metal Tall Fillers	11.04
Metal Knee Space, Panels & Rails.....	11.05
Metal Utility Tables	11.06-11.07

Base Fillers, Panels and Accessories



CAT. NO

R1200

DESCRIPTION

Base Scribing Filler

HEIGHT: INCHES

29 (Sitting) 33 (ADA) 35 (Standing)

WIDTH: INCHES

8 14



R1201

Base Front Filler

35 (Standing)



R1210

Base Chase Filler

29 (Sitting) 33 (ADA) 35 (Standing)

7



CAT. NO

R1211

DESCRIPTION


Base End Filler

- Used between back to back cabinets in an island configuration




R1250

Base Corner Filler



R2100

End Gable Panel



CAT. NO

PWS00

DESCRIPTION

Pull-Out Writing Shelf

- For sitting-height cabinet only

HEIGHT: INCHES


22

DEPTH: INCHES

15 18 21 24

WIDTH: INCHES

15 18 21 24




SHP00

Pull-Out Shelf for Base Cabinet

22

15 18 24 30 36 42 48

15 18 24 30 36 42 48



SHB00

Adjustable Shelf for Base Cabinet

22




15 18 24 30 36 42 48

15 18 24 30 36 42 48



CiFLAB
SOLUTIONS

Wall Fillers





CAT. NO	R1000	R1050	SHW00
DESCRIPTION	Wall Scribing Filler	Wall Corner Filler	Adjustable Shelf for Wall Case
HEIGHT: INCHES	24 31 36 48	24 31 36 48	
DEPTH: INCHES	13 16	13 16	22
WIDTH: INCHES	Specify	Specify	18 24 30 36 42 48





CiFLAB
SOLUTIONS

Tall Fillers




CAT. NO	R1300	R1350	SHT00
DESCRIPTION	Tall Cabinet Front Filler Pane	Tall Cabinet Corner Filler Panel With Metal Kick Plate • Shipped unassembled	Adjustable Shelf for Tall Case
HEIGHT: INCHES	84	84 (2134 MM)	
DEPTH: INCHES		12 15 22 (305 381 559 MM)	16 21
WIDTH: INCHES	2 8	1 2 6 12 (25 51 152 305 MM)	18 24 30 36 42 48

Knee Space, Panels and Rails



CAT. NO	R2201	R2200	R3300
DESCRIPTION	Knee Space Panel	Two Piece Knee Space Panel (Removable)	Knee Space Rail
HEIGHT: INCHES	24 (Sitting w/apron) 29 (Sitting without apron or standing w/apron) 35	29 35	4
WIDTH: INCHES	24 30 36 42 48 54 60 66 72	24 30 36 42 48 54 60 66 72	18 24 30 36 42 48 52 60 66 72 84 96

Metal Utility Tables



CAT. NO

D6140

DESCRIPTION

Table Apron

HEIGHT: INCHES


04

DEPTH: INCHES

22 28

WIDTH: INCHES

24 30 36 42 48 60 66 72




D6141

Table Apron with One Drawer

04

22 28

24 30 36 42




D6142

Table Apron with Two Drawers

04

22 28

42 48 54 60 66 72



CAT. NO

D6150

DESCRIPTION

Table Frame and Four Legs

HEIGHT: INCHES


29 (Sitting) 33 (ADA) 35 (Standing)

DEPTH: INCHES

22 28

WIDTH: INCHES

24 30 36 42 48 60 72




D6151

Table Frame, One Drawer and Four Legs

29 (Sitting) 33 (ADA) 35 (Standing)

22 28

24 30 36 42




D6152

Table Frame, Two Drawers and Four Legs

29 (Sitting) 33 (ADA) 35 (Standing)

22 28

42 48 54 60 66 72



CAT. NO

D6153

DESCRIPTION

Table Frame, Three Drawers and Four Legs

HEIGHT: INCHES

29 (Sitting) 35 (Standing)

DEPTH: INCHES

22 28

WIDTH: INCHES

72



CiFLAB
SOLUTIONS

Metal Utility Tables



CAT. NO	D6191	D6192	D6193
DESCRIPTION	Single Leg	Double Leg (Welded)	Single Leg w/Plate
HEIGHT: INCHES	29 (Sitting) 33 (ADA) 35 (Standing)	29 (Sitting) 33 (ADA) 35 (Standing)	29 (Sitting) 33 (ADA) 35 (Standing)
DEPTH: INCHES	22 28	22 28	
WIDTH: INCHES	2	2	



CiFLAB
SOLUTIONS

Contents

Metal Suspended & Mobile Base Cabinets	12.02
---	--------------



CiFLAB
SOLUTIONS

General Storage Mobile Cabinets



CAT. NO

U4030

DESCRIPTION

Storage With Three Drawers

- Does not include top

HEIGHT: INCHES

29

DEPTH: INCHES

22

WIDTH: INCHES

15 18 24 30

All metal base cabinets except sink, corner, acid, flammable and fume hood can be suspended or made mobile. Add the suffix as follows:

Suspended = HB

Mobile = MB