


# Contents


Rod Assemblies.....	AS-2
Computer Accessories.....	AS-3
Task Lights.....	AS-4
Miscellaneous.....	AS-5
Accessories Index .....	AS-6

# Rod Assemblies


**CAT. NO SA100**

**DESCRIPTION T-Bar Assembly**

- 1-SA050, 48" upright rod
- 1-SA060, 36" crossbar
- 1-SA010, adjustable clamp


**SA200**

**H-Bar Assembly**

- 2-SA050, 48" upright rods
- 1-SA060, 36" crossbar
- 2-SA010, adjustable clamps


**CAT. NO SA060**

**DESCRIPTION Crossbar**

- 36" length
- 3/4" aluminum crossbar


**SA070**

**Crossbar**

- 24" length
- 3/4" aluminum crossbar


**SA040**


**Upright Rod**

- 39" length
- 3/4" aluminum vertical rod with tapered end


**CAT. NO SA050**


**DESCRIPTION Upright Rod**

- 46" length
- 3/4" aluminum vertical rod with tapered end


**SA010**

**Adjustable Clamp**


- Aluminum 3/4" adjustable clamp










**SA020**

**Burette Socket**



- Aluminum burette socket with locknut
- Requires 1-1/8" hole in countertop

# Computer Accessories


**CAT. NO** **GR650**
**DESCRIPTION** Paper Grommet

**CAT. NO** **GR600**
**DESCRIPTION** Round Grommet

**CAT. NO** **GR620**
**DESCRIPTION** Oval Grommet

**CAT. NO** **RKB30**
**DESCRIPTION** Keyboard Tray w/ Mouse Tray

**CAT. NO** **RKB20**
**DESCRIPTION** Keyboard Tray

**CAT. NO** **CP030**
**DESCRIPTION** CPU Holder

**CAT. NO** **RKEY0**
**DESCRIPTION** Keyboard Tray

**CAT. NO** **RKB50**
**DESCRIPTION** Keyboard Tray w/ Mouse Tray

**CAT. NO** **WC010**
**DESCRIPTION** Wiring Tray

**CAT. NO** **WC020**
**DESCRIPTION** Wiring Tray

# Task Lights



CAT. NO **TL240**

**TLED0**


DESCRIPTION 24" Fluorescent Light

LED Task Light

# Miscellaneous




**CAT. NO** **DM010**  
**DESCRIPTION** **Demonstration Mirror**




**DC010**  
**Ceiling and Track Mounted Demonstration Mirror**




**LSSE1083**  
**First Aid Kit – 50 People**  
 • Designed to meet OSHA and ANSI first aid requirements.  
 • Durable, steel box has a built-in handle and is wall mountable.




**CAT. NO** **CC010**  
**DESCRIPTION** **Cylinder Clamp**



**VB040**  
**Black Vinyl Base**



**VC040**  
**Black Vinyl Base Corner**




**CAT. NO** **H8060**  
**DESCRIPTION** **Riser Shelf**  
 • Add to any work surface for additional storage  
 • Can be used for peripherals such as printers and scanners in order to free valuable work surface  
 • Available in oak or maple veneer or plastic laminate  
 • 12-1/2”h x 12”d x 24”w



**PC020**  
**4-sided (Island) Pipe Chase**  

<b>H</b>	48-96
<b>D</b>	03   06   09   12   15   18   21   24
<b>W</b>	04   06   08   10   12   14   16   18   20   22   24   26   28   30



**PC020**  
**3-Sided (Wall) Pipe Chase**  

<b>H</b>	48-96
<b>D</b>	03   06   09   12   15   18   21   24
<b>W</b>	04   06   08   10   12   14   16   18   20   22   24   26   28   30

# Accessories Index

<b>Accessories .....</b>	<b>AS-2</b>
SA100 T-Bar Assembly.....	AS-2
SA200 H-Bar Assembly .....	AS-2
SA060 Crossbar .....	AS-2
SA070 Crossbar .....	AS-2
SA040 Upright Rod.....	AS-2
SA050 Upright Rod.....	AS-2
SA010 Adjustable Clamp.....	AS-2
SA020 Burette Socket .....	AS-2
GR650 Paper Grommet .....	AS-3
GR600 Round Grommet .....	AS-3
GR620 Oval Grommet.....	AS-3
RKB30 Keyboard Tray w/ Mouse Tray .....	AS-3
RKB20 Keyboard Tray.....	AS-3
CP030 CPU Holder.....	AS-3
RKEY0 Keyboard Tray.....	AS-3
RKB50 Keyboard Tray w/ Mouse Tray .....	AS-3
WC010 Wiring Tray .....	AS-3
WC020 Wiring Tray.....	AS-3
TL240 24" Fluorescent Light.....	AS-4
TLED0 LED Task Light.....	AS-4
DM010 Demonstration Mirror .....	AS-5
DC010 Ceiling and Track Mounted Demonstration Mirror .....	AS-5
LSSE1083 First Aid Kit – 50 People.....	AS-5
CC010 Cylinder Clamp .....	AS-5
VB040 Black Vinyl Base .....	AS-5
VC040 Black Vinyl Base Corner .....	AS-5
H8060 Riser Shelf .....	AS-5
PC020 4-sided (Island) Pipe Chase.....	AS-5
PC020 3-Sided (Wall) Pipe Chase.....	AS-5