

**Project Name: Sharon City School District STEAM****Project Location: Sharon, PA****Project Overview**

Sharon City School District renovated the science classrooms within their Middle and High School facilities to accommodate STEAM curriculum. The renovation will provide students with a collaborative environment for learning and discovery. Northeast Interior Systems installed CiF Lab Solutions' Educator Series premium Tier 1 wood casework. The Educator Series incorporates integral end panels and full flush overlay 1/16" reveals. Red oak veneer with a clear finish and vertical grain match, for consistent uniformity in grain and color, complete the space.

- CiF Premium Educator Series Red Oak Wood Casework
- CiF Apex Air Fume Hoods
- Durcon Black Epoxy Countertops

Project Statistics

- Project Size: 18,568 sq. ft.
- Total Construction Budget: \$4.1 million
- CiF Contract: \$105,000

Project Participants

- Architect: Eckles Architecture
- General Contractor: Hudson Construction

