

# Contents

General Information.....	TB-2
Wood Tables Summary .....	TB-3
Fixed-Height Wood Tables With Leveling Glides.....	TB-4
Fixed-Height Wood Tables With Casters.....	TB-5
Adjustable-Height Wood Tables with Manual Adjustment and Leveling Glides .....	TB-6
Adjustable-Height Wood Tables with Manual Adjustment and Casters .....	TB-7
Wood Tables.....	TB-8
Metal Tables Summary .....	TB-9
Metal Tables .....	TB-10
Tables Index .....	TB-11

# General Information

Please review the following for general cabinet information. Contact your dealer for any clarifications.

## All wood tables

- Sizes shown are overall sizes including 1" thick top and 1" overhang all around.
- Solid hardwood (plain-sliced red oak or white maple) construction.
- Internal intermediate rails are plywood.
- Aprons are 4-5/8" high for standing-height, 3-3/4" high for ADA and sitting height.
- For tables over 72" wide, we recommend the use of one four-leg table with added left or right leg tables.
- We recommend a maximum length of 96" for multiple-leg tables for ease of movement and to allow making the top in one piece.

## Fixed-height wood tables with levelers

- Legs are 2-1/4" square hardwood and include black covered leg shoes.
- Tables over 48" wide include 1-1/4" x 1-3/4" leg stretchers for additional leg support consisting of one at each end and one longitudinal (specify center or rear mounted stretcher).

## Fixed-height wood tables with casters

- Legs are 2-1/4" square hardwood and include 100 lb. swivel casters (2 locking).
- All tables with casters include 1-1/4" x 1-3/4" leg stretchers for additional leg support consisting of one at each end and one longitudinal (specify center or rear mounted stretcher).

## Adjustable-height wood tables with levelers

- Legs are 2-1/2" square hardwood and include stainless steel telescoping inserts, adjustable from 30" to 37" high (by means of multiple holes & threaded thru-bolt) and leveler.
- Tables over 48" wide include 1-1/4" x 1-3/4" leg stretchers for additional leg support consisting of one at each end and one longitudinal (specify center or rear mounted stretcher).


## Adjustable-height wood tables with casters

- Legs are 2-1/2" square hardwood and include stainless steel telescoping inserts, adjustable from 30" to 37" high (by means of multiple holes & threaded thru-bolt) and 100 lb. swivel casters (2 locking).
- All tables with casters include 1-1/4" x 1-3/4" leg stretchers for additional leg support consisting of one at each end and one longitudinal (specify center or rear mounted stretcher).

# Wood Tables Summary


Description	No Top	Phenolic Top	Epoxy Resin Top	Maple Butcher Block Top
<b>Fixed-height with levelers</b>				
Table frame, 4 legs	D6150	D6150PR	D6150ER	D6150BB
Table frame, 1 drawer, 4 legs	D6151	D6151PR	D6151ER	D6151BB
Table frame, 2 drawers, 4 legs	D6152	D6152PR	D6152ER	D6152BB
Table frame, right legs	D6170	D6170PR	D6170ER	D6170BB
Table frame, 1 drawer, right legs	D6171	D6171PR	D6171ER	D6171BB
Table frame, 2 drawers, right legs	D6172	D6172PR	D6172ER	D6172BB
Table frame, left legs	D6160	D6160PR	D6160ER	D6160BB
Table frame, 1 drawer, left legs	D6161	D6161PR	D6161ER	D6161BB
Table frame, 2 drawers, left legs	D6162	D6162PR	D6162ER	D6162BB
Trapezoidal table	D6190	D6190PR	D6190ER	D6190BB
<b>Fixed-height with casters</b>				
Table frame, 4 legs, casters	D6250	D6250PR	D6250ER	D6250BB
Table frame, 1 drawer, 4 legs, casters	D6251	D6251PR	D6251ER	D6251BB
Table frame, 2 drawers, 4 legs, casters	D6252	D6252PR	D6252ER	D6252BB
Table frame, right legs, casters	D6270	D6270PR	D6270ER	D6270BB
Table frame, 1 drawer, right legs, casters	D6271	D6271PR	D6271ER	D6271BB
Table frame, 2 drawers, right legs, casters	D6272	D6272PR	D6272ER	D6272BB
Table frame, left legs, casters	D6260	D6260PR	D6260ER	D6260BB
Table frame, 1 drawer, left legs, casters	D6261	D6261PR	D6261ER	D6261BB
Table frame, 2 drawers, left legs, casters	D6262	D6262PR	D6262ER	D6262BB
Trapezoidal table, casters	D6290	D6290PR	D6290ER	D6290BB
<b>Adjustable-height (manual) with levelers</b>				
Table frame, 4 legs, adjustable height	D6350	D6350PR	D6350ER	D6350BB
Table frame, 1 drawer, 4 legs, adjustable height	D6351	D6351PR	D6351ER	D6351BB
Table frame, 2 drawers, 4 legs, adjustable height	D6352	D6352PR	D6352ER	D6352BB
Table frame, right legs, adjustable height	D6370	D6370PR	D6370ER	D6370BB
Table frame, 1 drawer, right legs, adjustable height	D6371	D6371PR	D6371ER	D6371BB
Table frame, 2 drawers, right legs, adjustable height	D6372	D6372PR	D6372ER	D6372BB
Table frame, left legs, adjustable height	D6360	D6360PR	D6360ER	D6360BB
Table frame, 1 drawer, left legs, adjustable height	D6361	D6361PR	D6361ER	D6361BB
Table frame, 2 drawers, left legs, adjustable height	D6362	D6362PR	D6362ER	D6362BB
Trapezoidal table, adjustable height	D6390	D6390PR	D6390ER	D6390BB
<b>Adjustable-height (manual) with casters</b>				
Table frame, 4 legs, adjustable height, casters	D6450	D6450PR	D6450ER	D6450BB
Table frame, 1 drawer, 4 legs, adjustable height, casters	D6451	D6451PR	D6451ER	D6451BB
Table frame, 2 drawers, 4 legs, adjustable height, casters	D6452	D6452PR	D6452ER	D6452BB
Table frame, right legs, adjustable height, casters	D6470	D6470PR	D6470ER	D6470BB
Table frame, 1 drawer, right legs, adjustable height, casters	D6471	D6471PR	D6471ER	D6471BB
Table frame, 2 drawers, right legs, adjustable height, casters	D6472	D6472PR	D6472ER	D6472BB
Table frame, left legs, adjustable height, casters	D6460	D6460PR	D6460ER	D6460BB
Table frame, 1 drawer, left legs, adjustable height, casters	D6461	D6461PR	D6461ER	D6461BB
Table frame, 2 drawers, left legs, adjustable height, casters	D6462	D6462PR	D6462ER	D6462BB
Trapezoidal table, adjustable height, casters	D6490	D6490PR	D6490ER	D6490BB
<b>Apron frames only (to be used in conjunction with fixed casework)</b>				
Table frame, no legs	D6140			
Table frame, 1 drawer, no legs	D6141			
Table frame, 2 drawers, no legs	D6142			
Table leg only, end	D6191			

# Fixed-Height Wood Tables With Leveling Glides


**CAT. NO** **D6150**
**DESCRIPTION** Table Frame


**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6151**
**DESCRIPTION** Table Frame with One Drawer


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6152**
**DESCRIPTION** Table Frame with Two Drawers


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6170**
**DESCRIPTION** Table Frame – Right Legs


**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6171**
**DESCRIPTION** Table Frame with One Drawer – Right Legs


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6172**
**DESCRIPTION** Table Frame with Two Drawers – Right Legs


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6160**
**DESCRIPTION** Table Frame – Left Legs


**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6161**
**DESCRIPTION** Table Frame with One Drawer – Left Legs


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6162**
**DESCRIPTION** Table Frame with Two Drawers – Left Legs

**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30


**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6190**
**DESCRIPTION** Trapezoidal Table Frame

**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30


**WIDTH** 30 | 60

# Fixed-Height Wood Tables With Casters


**CAT. NO** **D6250**
**DESCRIPTION** Table Frame


**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6251**
**DESCRIPTION** Table Frame with One Drawer


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6252**
**DESCRIPTION** Table Frame with Two Drawers


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6270**
**DESCRIPTION** Table Frame – Right Legs


**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6271**
**DESCRIPTION** Table Frame with One Drawer – Right Legs


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6272**
**DESCRIPTION** Table Frame with Two Drawers – Right Legs


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6260**
**DESCRIPTION** Table Frame – Left Legs


**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6261**
**DESCRIPTION** Table Frame with One Drawer – Left Legs


**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6262**
**DESCRIPTION** Table Frame with Two Drawers – Left Legs

**HEIGHT** 30 (sitting) | 36 (standing)

**DEPTH** 24 | 30


**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6290**
**DESCRIPTION** Trapezoidal Table Frame

**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30


**WIDTH** 30 | 60

# Adjustable-Height Wood Tables with Manual Adjustment and Leveling Glides


**CAT. NO** **D6350**
**DESCRIPTION** Table Frame


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6351**
**DESCRIPTION** Table Frame with One Drawer


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6352**
**DESCRIPTION** Table Frame with Two Drawers


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6370**
**DESCRIPTION** Table Frame – Right Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6371**
**DESCRIPTION** Table Frame with One Drawer – Right Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6372**
**DESCRIPTION** Table Frame with Two Drawers – Right Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6360**
**DESCRIPTION** Table Frame – Left Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6361**
**DESCRIPTION** Table Frame with One Drawer – Left Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6362**
**DESCRIPTION** Table Frame with Two Drawers – Left Legs

**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30


**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6390**
**DESCRIPTION** Trapezoidal Table Frame

**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30


**WIDTH** 30 | 60

# Adjustable-Height Wood Tables with Manual Adjustment and Casters


**CAT. NO** **D6450**
**DESCRIPTION** Table Frame


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6451**
**DESCRIPTION** Table Frame with One Drawer


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6452**
**DESCRIPTION** Table Frame with Two Drawers


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6470**
**DESCRIPTION** Table Frame – Right Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6471**
**DESCRIPTION** Table Frame with One Drawer – Right Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6472**
**DESCRIPTION** Table Frame with Two Drawers – Right Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6460**
**DESCRIPTION** Table Frame – Left Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30 | 48

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **D6461**
**DESCRIPTION** Table Frame with One Drawer – Left Legs


**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36

**CAT. NO** **D6462**
**DESCRIPTION** Table Frame with Two Drawers – Left Legs

**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 42 | 48 | 60 | 72

**CAT. NO** **D6490**
**DESCRIPTION** Trapezoidal Table Frame




**HEIGHT** 30 - 37" Adjustable

**DEPTH** 24 | 30

**WIDTH** 30 | 60

# Wood Tables


Table frames on this page are intended for use in conjunction with fixed casework. Tops are not included.



<b>CAT. NO</b>	<b>D6140</b>
<b>DESCRIPTION</b>	Table Frame – No Legs
<b>HEIGHT – in</b>	3-3/4 (Sitting & ADA)   4-5/8 (Standing)
<b>HEIGHT – mm</b>	95   117
<b>DEPTH – in</b>	22   28
<b>DEPTH – mm</b>	559   711
<b>WIDTH – in</b>	24   30   36   42   48   60   72
<b>WIDTH – mm</b>	610   762   914   1067   1219   1524   1829

<b>CAT. NO</b>	<b>D6141</b>
<b>DESCRIPTION</b>	Table Frame with One Drawer – No Legs
<b>HEIGHT – in</b>	3-3/4 (Sitting & ADA)   4-5/8 (Standing)
<b>HEIGHT – mm</b>	95   117
<b>DEPTH – in</b>	22   28
<b>DEPTH – mm</b>	559   711
<b>WIDTH – in</b>	24   30   36
<b>WIDTH – mm</b>	610   762   914

<b>CAT. NO</b>	<b>D6142</b>
<b>DESCRIPTION</b>	Table Frame with Two Drawers – No Legs
<b>HEIGHT – in</b>	3-3/4 (Sitting & ADA)   4-5/8 (Standing)
<b>HEIGHT – mm</b>	95   117
<b>DEPTH – in</b>	22   28
<b>DEPTH – mm</b>	559   711
<b>WIDTH – in</b>	42   48   60   72
<b>WIDTH – mm</b>	1067   1219   1524   1829



<b>CAT. NO</b>	<b>D6191</b>
<b>DESCRIPTION</b>	Table Leg – End
<b>HEIGHT – in</b>	29 (sitting)   32-3/4 (ADA)   35 (standing)
<b>HEIGHT – mm</b>	737   832   889



# General Information – Metal Tables

Please review the following for general information. Contact your dealer for any clarifications.

## All metal tables

- Sizes shown are overall sizes including 1" thick top and 1" overhang all around.
- Frames and legs are 2" x 2" 16 gauge steel with all-welded construction.
- Epoxy powder-coated finish (black is standard).

## Fixed-height metal tables with levelers

- Equipped with levelers and metal leg shoes.
- Tables over 48" wide include leg stretchers for additional leg support consisting of one at each end and one longitudinal (specify center or rear mounted stretcher).

## Fixed-height metal tables with casters

- Include 100 lb. swivel casters (2 locking).
- All tables with casters include leg stretchers for additional leg support consisting of one at each end and one longitudinal (specify center or rear mounted stretcher).

## Adjustable-height metal tables

- Adjustable in 1" increments between 26" high and 38" high. Legs include painted metal inserts which are manually adjustable by means of a bolt through the leg to fix the height. Includes levelers and metal leg shoes.
- All tables with adjustable-height include leg stretchers for additional leg support consisting of one at each end and one longitudinal (specify center or rear mounted stretcher).

## Adjustable-height metal tables with casters

- Adjustable in 1" increments between 26" high and 38" high. Legs include painted metal inserts which are manually adjustable by means of a bolt through the leg to fix the height. Includes levelers and metal leg shoes.
- Include 100 lb. swivel casters (2 locking).
- All tables with adjustable-height include leg stretchers for additional leg support consisting of one at each end and one longitudinal (specify center or rear mounted stretcher).


Description	No Top	Phenolic Top	Epoxy Resin Top	Maple Butcher Block
Metal table, fixed height with levelers	V6150	V6150PR	V6150ER	V6150BB
Metal table, fixed height with casters	V6250	V6250PR	V6250ER	V6250BB
Metal table, adjustable-height with levelers	V6350	V6350PR	V6350ER	V6350BB
Metal table, adjustable-height with casters	V6450	V6450PR	V6450ER	V6450BB
Metal trapezoidal table, fixed height with levelers	V6190	V6190PR	V6190ER	V6190BB
Metal trapezoidal table, fixed height with casters	V6290	V6290PR	V6290ER	V6290BB
Metal trapezoidal table, adjustable-height with levelers	V6390	V6390PR	V6390ER	V6390BB
Metal trapezoidal table, adjustable-height with casters	V6490	V6490PR	V6490ER	V6490BB

# Metal Tables


**CAT. NO** **V6150**
**DESCRIPTION** Metal table, fixed height with levelers


**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **V6190**
**DESCRIPTION** Metal trapezoidal table, fixed height with levelers


**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 48 | 60

**CAT. NO** **V6250**
**DESCRIPTION** Metal table, fixed height with casters

**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **V6290**
**DESCRIPTION** Metal trapezoidal table, fixed height with casters

**HEIGHT** 30 (sitting) | 34 (ADA) | 36 (standing)


**DEPTH** 24 | 30

**WIDTH** 48 | 60

**CAT. NO** **V6350**
**DESCRIPTION** Metal table, adjustable-height with levelers


**HEIGHT** 26 - 38" Adjustable

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **V6390**
**DESCRIPTION** Metal trapezoidal table, adjustable-height with levelers


**HEIGHT** 26 - 38" Adjustable

**DEPTH** 24 | 30

**WIDTH** 48 | 60

**CAT. NO** **V6450**
**DESCRIPTION** Metal table, adjustable-height with casters

**HEIGHT** 26 - 38" Adjustable

**DEPTH** 24 | 30

**WIDTH** 24 | 30 | 36 | 42 | 48 | 60 | 72

**CAT. NO** **V6490**
**DESCRIPTION** Metal trapezoidal table, adjustable-height with casters

**HEIGHT** 26 - 38" Adjustable

**DEPTH** 24 | 30

**WIDTH** 48 | 60

# Tables Index

<b>General Information</b> .....	<b>TB-2</b>	D6360	Table Frame – Left Legs .....	TB-6	
<b>Wood Tables Summary</b> .....	<b>TB-3</b>	D6361	Table Frame with One Drawer – Left Legs .....	TB-6	
<b>Fixed-Height Wood Tables With Leveling Glides</b> .....	<b>TB-4</b>	D6362	Table Frame with Two Drawers – Left Legs .....	TB-6	
D6150	Table Frame .....	TB-4	D6390	Trapezoidal Table Frame .....	TB-6
D6151	Table Frame with One Drawer .....	TB-4	<b>Adjustable-Height Wood Tables with Manual Adjustment and Casters</b> .....		
D6152	Table Frame with Two Drawers .....	TB-4	<b>TB-7</b>		
D6170	Table Frame – Right Legs .....	TB-4	D6450	Table Frame .....	TB-7
D6171	Table Frame with One Drawer – Right Legs .....	TB-4	D6451	Table Frame with One Drawer .....	TB-7
D6172	Table Frame with Two Drawers – Right Legs .....	TB-4	D6452	Table Frame with Two Drawers .....	TB-7
D6160	Table Frame – Left Legs .....	TB-4	D6470	Table Frame – Right Legs .....	TB-7
D6161	Table Frame with One Drawer – Left Legs .....	TB-4	D6471	Table Frame with One Drawer – Right Legs .....	TB-7
D6162	Table Frame with Two Drawers – Left Legs .....	TB-4	D6472	Table Frame with Two Drawers – Right Legs .....	TB-7
D6190	Trapezoidal Table Frame .....	TB-4	D6460	Table Frame – Left Legs .....	TB-7
<b>Fixed-Height Wood Tables With Casters</b> .....	<b>TB-5</b>	D6461	Table Frame with One Drawer – Left Legs .....	TB-7	
D6250	Table Frame .....	TB-5	D6462	Table Frame with Two Drawers – Left Legs .....	TB-7
D6251	Table Frame with One Drawer .....	TB-5	D6490	Trapezoidal Table Frame .....	TB-7
D6252	Table Frame with Two Drawers .....	TB-5	<b>Wood Tables</b> .....		
D6270	Table Frame – Right Legs .....	TB-5	<b>TB-8</b>		
D6271	Table Frame with One Drawer – Right Legs .....	TB-5	D6140	Table Frame – No Legs .....	TB-8
D6272	Table Frame with Two Drawers – Right Legs .....	TB-5	D6141	Table Frame with One Drawer – No Legs .....	TB-8
D6260	Table Frame – Left Legs .....	TB-5	D6142	Table Frame with Two Drawers – No Legs .....	TB-8
D6261	Table Frame with One Drawer – Left Legs .....	TB-5	D6191	Table Leg – End .....	TB-8
D6262	Table Frame with Two Drawers – Left Legs .....	TB-5	<b>Metal Tables Summary</b> .....		
D6290	Trapezoidal Table Frame .....	TB-5	<b>TB-9</b>		
<b>Adjustable-Height Wood Tables with Manual Adjustment and Leveling Glides</b> .....	<b>TB-6</b>	<b>Metal Tables</b> .....		<b>TB-10</b>	
D6350	Table Frame .....	TB-6	V6150	Metal table, fixed height with levelers .....	10
D6351	Table Frame with One Drawer .....	TB-6	V6190	Metal trapezoidal table, fixed height with levelers .....	10
D6352	Table Frame with Two Drawers .....	TB-6	V6250	Metal table, fixed height with casters .....	10
D6370	Table Frame – Right Legs .....	TB-6	V6290	Metal trapezoidal table, fixed height with casters .....	10
D6371	Table Frame with One Drawer – Right Legs .....	TB-6	V6350	Metal table, adjustable-height with levelers .....	10
D6372	Table Frame with Two Drawers – Right Legs .....	TB-6	V6390	Metal trapezoidal table, adjustable-height with levelers .....	10
			V6450	Metal table, adjustable-height with casters .....	10
			V6490	Metal trapezoidal table, adjustable-height with casters .....	10

